

ELOVL4 protein (STJ504004)

STJ504004

GENERAL INFORMATION

Product Type	Proteins
Short Description	ELOVL4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6785
Gene Symbol	ELOVL4
Uniprot ID	ELOV4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 294-301 of Human ELOVL4. Sequence: SEKQLMIENGKQKNGKAKGDC
Immunogen Region	294-301
Specificity	
Immunogen Sequence	