

Prostaglandin E Receptor EP2 protein (STJ504007)

STJ504007

GENERAL INFORMATION

Product Type	Proteins
Short Description	Prostaglandin E Receptor EP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5733
Gene Symbol	PTGER3
Uniprot ID	PE2R3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Prostaglandin E2 receptor EP3 subtype.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	