

EPAC1 protein (STJ504008)

STJ504008

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | EPAC1-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 10411 |
| Gene Symbol | RAPGEF3 |
| Uniprot ID | RPGF3_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 5-50 on Exchange Protein Directly Activated by cAMP protein. |
| Immunogen Region | 5-50 |
| Specificity | |
| Immunogen Sequence | |