

## Ephrin Receptor A4 protein (STJ504012)

STJ504012

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Ephrin Receptor A4-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	13838
<b>Gene Symbol</b>	Epha4
<b>Uniprot ID</b>	EPHA4_MOUSE
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 550-650 on human EphA R4 protein.
<b>Immunogen Region</b>	550-650
<b>Specificity</b>	
<b>Immunogen Sequence</b>	