

## Ephrin Receptor B3 protein (STJ504019)

STJ504019

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Ephrin Receptor B3-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">13845</a>
<b>Gene Symbol</b>	<a href="#">Ephb3</a>
<b>Uniprot ID</b>	<a href="#">EPHB3_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 100-200 on EphBR3 protein.
<b>Immunogen Region</b>	100-200
<b>Specificity</b>	
<b>Immunogen Sequence</b>	