

Ephrin Receptor B4 protein (STJ504020)

STJ504020

GENERAL INFORMATION

Product Type	Proteins
Short Description	Ephrin Receptor B4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2050
Gene Symbol	EPHB4
Uniprot ID	EPHB4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-987 on human Ephrin B4 receptor protein.
Immunogen Region	900-987
Specificity	
Immunogen Sequence	