

ERG protein (STJ504021)

STJ504021

GENERAL INFORMATION

Product Type	Proteins
Short Description	ERG-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2078
Gene Symbol	ERG
Uniprot ID	ERG_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Transcriptional regulator ERG protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	