

HE4 protein (STJ504026)

STJ504026

GENERAL INFORMATION

Product Type	Proteins
Short Description	HE4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10406
Gene Symbol	WFDC2
Uniprot ID	WFDC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 70-120 on human WAP four-disulfide core domain protein 2.
Immunogen Region	70-120
Specificity	
Immunogen Sequence	