

Estrogen Receptor alpha protein (STJ504027)

STJ504027

GENERAL INFORMATION

Product Type	Proteins
Short Description	Estrogen Receptor alpha-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2099
Gene Symbol	ESR1
Uniprot ID	ESR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 30-80 on human ESTRA protein
Immunogen Region	30-80
Specificity	
Immunogen Sequence	