

ergo-1 protein (STJ504058)
STJ504058**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ergo-1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 790-840 on C.elegans Endogenous-RNAi deficient arGOnaute protein.
Immunogen Region	790-840
Specificity	
Immunogen Sequence	