

FAK1 protein (STJ504061)

STJ504061

GENERAL INFORMATION

Product Type	Proteins
Short Description	FAK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5747
Gene Symbol	PTK2
Uniprot ID	FAK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 740-790 on human focal adhesion kinase 1 protein.
Immunogen Region	740-790
Specificity	
Immunogen Sequence	