

FAM38A protein (STJ504063)

STJ504063

GENERAL INFORMATION

Product Type	Proteins
Short Description	FAM38A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9780
Gene Symbol	PIEZO1
Uniprot ID	PIEZ1_HUMAN
Immunogen	Synthetic peptide taken between TMD14 and TMD15 on human FAM38A protein.
Immunogen Region	
Specificity	
Immunogen Sequence	