

FBXW4 protein (STJ504064)
STJ504064**GENERAL INFORMATION**

Product Type	Proteins
Short Description	FBXW4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6468
Gene Symbol	FBXW4
Uniprot ID	FBXW4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 240-290 on human FBX4 protein.
Immunogen Region	240-290
Specificity	
Immunogen Sequence	