

GPCR GPR43 protein (STJ504070)

STJ504070

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR43-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2867
Gene Symbol	FFAR2
Uniprot ID	FFAR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-300 on human free fatty acid receptor 2 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	