

GPCR GPR120 protein (STJ504073)

STJ504073

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR120-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	338557
Gene Symbol	FFAR4
Uniprot ID	FFAR4_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-80 on human FFAR4 protein.
Immunogen Region	1-80
Specificity	
Immunogen Sequence	