

GPCR GPR120 FFAR4-412AP protein (STJ504074)

STJ504074

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR120 FFAR4-412AP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	338557
Gene Symbol	FFAR4
Uniprot ID	FFAR4_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 220-300 on human FFAR4 protein.
Immunogen Region	220-300
Specificity	
Immunogen Sequence	