

FLI1 protein (STJ504080)
STJ504080**GENERAL INFORMATION**

Product Type	Proteins
Short Description	FLI1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2313
Gene Symbol	FLI1
Uniprot ID	FLI1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Friend Leukemia integration 1 transcription factor protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	