

C11orf1 FLJ-112AP protein (STJ504083)

STJ504083

GENERAL INFORMATION

Product Type	Proteins
Short Description	C11orf1 FLJ-112AP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64776
Gene Symbol	C11orf1
Uniprot ID	CK001_HUMAN
Immunogen	Mid-region epitope synthetic peptide taken within amino acid region 50-100 on human C11orf1 protein
Immunogen Region	50-100
Specificity	
Immunogen Sequence	