

FOXP3 protein (STJ504085)

STJ504085

GENERAL INFORMATION

Product Type	Proteins
Short Description	FOXP3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64776
Gene Symbol	C11orf1
Uniprot ID	CK001_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 229-250 of Mouse Forkhead Box Protein 3 (winged helix transcription factor 3 protein). The peptide was post synthetically modified to achieve desired antigenicity before coupling to a carrier protein. See
Immunogen	229-250
Region	
Specificity	
Immunogen	
Sequence	