

FUBP1 protein (STJ504088)

STJ504088

GENERAL INFORMATION

Product Type	Proteins
Short Description	FUBP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8880
Gene Symbol	FUBP1
Uniprot ID	FUBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-644 on human FUBP1 protein.
Immunogen Region	600-644
Specificity	
Immunogen Sequence	