

FXIII A protein (STJ504091)

STJ504091

GENERAL INFORMATION

Product Type	Proteins
Short Description	FXIII A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2162
Gene Symbol	F13A1
Uniprot ID	F13A_HUMAN
Immunogen	Multi-epitope synthetic peptides were used: amino acid range 1-50 and 700-732. A mixture of synthetic peptides used for FXIII A-101AP and FXIII A-112AP.
Immunogen Region	1-50/700-732
Specificity	
Immunogen Sequence	