

G0/G1 switch 2 protein (STJ504116)

STJ504116

GENERAL INFORMATION

Product Type	Proteins
Short Description	G0/G1 switch 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 78-103 on mouse G0S2-like protein. The peptide was covalently modified to achieve desired antigenicity and conjugated to a carrier protein. Sequence: ALLAGGKAQEATLCSRALSLRQHAS
Immunogen Region	78-103
Specificity	
Immunogen Sequence	