

GABA B Receptor 2 protein (STJ504118)

STJ504118

GENERAL INFORMATION

Product Type	Proteins
Short Description	GABA B Receptor 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2561
Gene Symbol	GABRB2
Uniprot ID	GBRB2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human gamma-aminobutyric acid A receptor beta 2 protein. Will label various GABA B subunit 2 variants.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	