

GABA B Receptor 2 protein (STJ504119)

STJ504119

GENERAL INFORMATION

Product Type	Proteins
Short Description	GABA B Receptor 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2561
Gene Symbol	GABRB2
Uniprot ID	GBRB2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human gamma-aminobutyric acid receptor subunit beta-2 protein.
Immunogen Region	350-400
Specificity	
Immunogen	
Sequence	