

GADD34 protein (STJ504123)

STJ504123

GENERAL INFORMATION

Product Type	Proteins
Short Description	GADD34-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to internal sequence of amino acids 380-400 of Human GADD34, covalently modified and coupled to a carrier protein. Sequence: PATGVFLKSWVYQPGEDTEEE
Immunogen Region	380-400
Specificity	
Immunogen Sequence	