

GCN2 protein (STJ504124)

STJ504124

GENERAL INFORMATION

Product Type	Proteins
Short Description	GCN2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	440275
Gene Symbol	EIF2AK4
Uniprot ID	E2AK4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1600-1700 on human GCN2 protein
Immunogen Region	1600-1700
Specificity	
Immunogen Sequence	