

GGCX protein (STJ504126)

STJ504126

GENERAL INFORMATION

Product Type	Proteins
Short Description	GGCX-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2677
Gene Symbol	GGCX
Uniprot ID	VKGC_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human GGCX protein
Immunogen Region	400-450
Specificity	
Immunogen Sequence	