

## Growth Hormone protein (STJ504128)

STJ504128

### GENERAL INFORMATION

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| <b>Product Type</b>      | Proteins                                                    |
| <b>Short Description</b> | Growth Hormone-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                          |
| <b>Host/Source</b>       |                                                             |
| <b>Reactivity</b>        |                                                             |

### PRODUCT PROPERTIES

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Clonality</b>           |                                                                 |
| <b>Clone ID</b>            |                                                                 |
| <b>Concentration</b>       |                                                                 |
| <b>Conjugation</b>         |                                                                 |
| <b>Purification</b>        |                                                                 |
| <b>Dilution Range</b>      |                                                                 |
| <b>Formulation</b>         |                                                                 |
| <b>Isotype</b>             |                                                                 |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            |                                                                                                                  |
| <b>Gene Symbol</b>        |                                                                                                                  |
| <b>Uniprot ID</b>         |                                                                                                                  |
| <b>Immunogen</b>          | C-epitope synthetic peptide taken within amino acid region 150-250 on human Growth Hormone beta subunit protein. |
| <b>Immunogen Region</b>   | 150-250                                                                                                          |
| <b>Specificity</b>        |                                                                                                                  |
| <b>Immunogen Sequence</b> |                                                                                                                  |