

Glutamate Receptor 1 protein (STJ504134)

STJ504134

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glutamate Receptor 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2890
Gene Symbol	GRIA1
Uniprot ID	GRIA1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-950.
Immunogen Region	850-950
Specificity	
Immunogen Sequence	