

Glutamate Receptor 2 protein (STJ504135)

STJ504135

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glutamate Receptor 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2891
Gene Symbol	GRIA2
Uniprot ID	GRIA2_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid region on GluR2 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	