

Glucose Transporter GLUT4 protein (STJ504144)

STJ504144

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glucose Transporter GLUT4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20528
Gene Symbol	Slc2a4
Uniprot ID	GLUT4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 450-509 on mouse Solute carrier family 2, facilitated glucose transporter member 4 protein.
Immunogen Region	450-509
Specificity	
Immunogen	
Sequence	