

Glucose transporter GLUT5 protein (STJ504146)

STJ504146

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glucose transporter GLUT5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6518
Gene Symbol	SLC2A5
Uniprot ID	GTR5_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human GLUT 5 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	