

Glucose Transporter 8 protein (STJ504148)

STJ504148

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glucose Transporter 8-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29988
Gene Symbol	SLC2A8
Uniprot ID	GTR8_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-500 on GLUT8 protein.
Immunogen Region	400-500
Specificity	
Immunogen Sequence	