

Golgin A1 protein (STJ504150)

STJ504150

GENERAL INFORMATION

Product Type	Proteins
Short Description	Golgin A1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	76899
Gene Symbol	Golga1
Uniprot ID	GOGA1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 650-700 on mouse Golgin subfamily A member 1 isoform b protein.
Immunogen Region	650-700
Specificity	
Immunogen Sequence	