

Golgin A2 protein (STJ504151)

STJ504151

GENERAL INFORMATION

Product Type	Proteins
Short Description	Golgin A2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 925-975 on human Golgin subfamily A member 2 isoform X3 protein.
Immunogen Region	925-975
Specificity	
Immunogen Sequence	