

Golgin A3 protein (STJ504152)

STJ504152

GENERAL INFORMATION

Product Type	Proteins
Short Description	Golgin A3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2802
Gene Symbol	GOLGA3
Uniprot ID	GOGA3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1400-1450 on human Golgin subfamily A member 3 isoform 1 protein.
Immunogen Region	1400-1450
Specificity	
Immunogen Sequence	