

Golgin A5 protein (STJ504154)

STJ504154

GENERAL INFORMATION

Product Type	Proteins
Short Description	Golgin A5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9950
Gene Symbol	GOLGA5
Uniprot ID	GOGA5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 660-710 on human Golgin subfamily A member 5 protein.
Immunogen Region	660-710
Specificity	
Immunogen Sequence	