

GPCR GPR3 protein (STJ504156)

STJ504156

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	266769
Gene Symbol	Gpr3
Uniprot ID	GPR3_RAT
Immunogen	Synthetic peptide taken within amino acid region 5-50 on rat GPCR3 protein.
Immunogen Region	5-50
Specificity	
Immunogen Sequence	