

GPCR GPR75 protein (STJ504159)

STJ504159

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR75-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10936
Gene Symbol	GPR75
Uniprot ID	GPR75_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 25-75 on human G protein-coupled receptor 75 protein.
Immunogen Region	25-75
Specificity	
Immunogen Sequence	