

GPCR GPR87 protein (STJ504161)

STJ504161

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR87-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53836
Gene Symbol	GPR87
Uniprot ID	GPR87_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on human GPCR87 protein.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	