

Granzyme B protein (STJ504165)

STJ504165

GENERAL INFORMATION

Product Type	Proteins
Short Description	Granzyme B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3002
Gene Symbol	GZMB
Uniprot ID	GRAB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human endogenous Granzyme B protein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	