

## GSK3 beta protein (STJ504166)

STJ504166

### GENERAL INFORMATION

|                          |                                                        |
|--------------------------|--------------------------------------------------------|
| <b>Product Type</b>      | Proteins                                               |
| <b>Short Description</b> | GSK3 beta-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                     |
| <b>Host/Source</b>       |                                                        |
| <b>Reactivity</b>        |                                                        |

### PRODUCT PROPERTIES

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Clonality</b>           |                                                                 |
| <b>Clone ID</b>            |                                                                 |
| <b>Concentration</b>       |                                                                 |
| <b>Conjugation</b>         |                                                                 |
| <b>Purification</b>        |                                                                 |
| <b>Dilution Range</b>      | WB:<br>IHC: 1:100                                               |
| <b>Formulation</b>         |                                                                 |
| <b>Isotype</b>             |                                                                 |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                                                                   |
|---------------------------|-----------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 2932                                                                              |
| <b>Gene Symbol</b>        | GSK3B                                                                             |
| <b>Uniprot ID</b>         | GSK3B_HUMAN                                                                       |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 1-50 on human GSK3 beta protein. |
| <b>Immunogen Region</b>   | 1-50                                                                              |
| <b>Specificity</b>        |                                                                                   |
| <b>Immunogen Sequence</b> |                                                                                   |