

Glucose-6-Phosphate, catalytic subunit protein (STJ504172)

STJ504172

GENERAL INFORMATION

Product Type	Proteins
Short Description	Glucose-6-Phosphate, catalytic subunit-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14378
Gene Symbol	G6pc2
Uniprot ID	G6PC2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	