

KCNJ2 protein (STJ504184)

STJ504184

GENERAL INFORMATION

Product Type	Proteins
Short Description	KCNJ2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:250
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3765
Gene Symbol	KCNJ9
Uniprot ID	KCNJ9_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	