

GPCR GPR3 protein (STJ504206)

STJ504206

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR GPR3-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:750
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	266769
Gene Symbol	Gpr3
Uniprot ID	GPR3_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	