

HAND2 protein (STJ504216)

STJ504216

GENERAL INFORMATION

Product Type	Proteins
Short Description	HAND2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64637
Gene Symbol	Hand2
Uniprot ID	HAND2_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on Heart and neuronal crest derivative protein 2.
Immunogen Region	
Specificity	
Immunogen Sequence	