

Helicobacter pylori protein (STJ504223)

STJ504223

GENERAL INFORMATION

Product Type	Proteins
Short Description	Helicobacter pylori-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10133
Gene Symbol	OPTN
Uniprot ID	OPTN_HUMAN
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Specificity	
Immunogen Sequence	