

ErbB 2 protein (STJ504226)
STJ504226**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ErbB 2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2064
Gene Symbol	ERBB2
Uniprot ID	ERBB2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	