

## Hepatitis C virus Core NS3 + NS4 protein (STJ504232)

STJ504232

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | Hepatitis C virus Core NS3 + NS4-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition  |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            |   |
| <b>Gene Symbol</b>        |   |
| <b>Uniprot ID</b>         |   |
| <b>Immunogen</b>          | N-epitope specific synthetic peptide corresponding to unique amino acid sequences on hepatitis C polyprotein. |
| <b>Immunogen Region</b>   | N-Term  |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |