

Hepatitis C virus Core NS3 + NS4 protein (STJ504233)

STJ504233

GENERAL INFORMATION

Product Type	Proteins
Short Description	Hepatitis C virus Core NS3 + NS4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	C-epitope specific synthetic peptide corresponding to unique amino acid sequences on hepatitis C polyprotein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	